

LAD1 antibody

Product Information

Catalog No.:	FNab04675
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Anchoring filament protein which is a component of the basement membrane zone.

Immunogen information


Immunogen:	ladinin 1
Synonyms:	Ladinin-1 (Lad-1) Linear IgA disease antigen (LADA) LAD1 LAD
Observed MW:	62 kDa
Uniprot ID :	O00515

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab04675(LAD1 antibody) at dilution of 1:100



human testis tissue were subjected to SDS PAGE followed by western blot with FNab04675(LAD1 antibody) at dilution of 1:3000