

## LAG3 antibody

### Product Information

Catalog No.:	FNab04676
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	lymphocyte-activation gene 3
Synonyms:	Lymphocyte activation gene 3 protein (LAG-3) Secreted lymphocyte activation gene 3 protein (sLAG-3) LAG3 FDC
Observed MW:	54 kDa, 60-65 kDa
Uniprot ID :	P18627

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	



IP Result of anti-LAG-3 (IP:FNab04676, 4ug;  
Detection:FNab04676 1:500) with mouse liver  
tissue lysate 4000ug.



mouse liver tissue were subjected to SDS PAGE  
followed by western blot with FNab04676(LAG-3  
antibody) at dilution of 1:500