

LAMP3 antibody

Product Information

Catalog No.:	FNab04687
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in dendritic cell function and in adaptive immunity.

Immunogen information


Immunogen:	lysosomal-associated membrane protein 3
Synonyms:	Lysosome-associated membrane glycoprotein 3 (LAMP-3, Lysosomal-associated membrane protein 3) DC-lysosome-associated membrane glycoprotein (DC LAMP) Protein TSC403 LAMP3 DCLAMP TSC403
Observed MW:	70 kDa
Uniprot ID :	Q9UQV4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNaB04687(LAMP3 antibody) at dilution of 1:50



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNaB04687(LAMP3 antibody) at dilution of 1:1000