

RPSA antibody

Product Information

Catalog No.:	FNab04688
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	ribosomal protein SA
Synonyms:	Small ribosomal subunit protein uS2 37 kDa laminin receptor precursor (37LRP) 37/67 kDa laminin receptor (LRP/LR) 40S ribosomal protein SA 67 kDa laminin receptor (67LR) Colon carcinoma laminin-binding protein Laminin receptor 1 (LamR) Laminin-binding protein precursor p40 (LBP/p40) Multidrug resistance-associated protein MGr1-Ag NEM/1CHD4 RPSA LAMBR LAMR1
Observed MW:	None
Uniprot ID :	P08865

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human ovary tumor using FNab04688(LAMR1,RPSA antibody) at dilution of 1:50