

LANCL1 antibody

Product Information

Catalog No.:	FNab04691
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in EPS8 signaling. Binds glutathione.

Immunogen information


Immunogen:	LanC lantibiotic synthetase component C-like 1(bacterial)
Synonyms:	Glutathione S-transferase LANCL1 40 kDa erythrocyte membrane protein (p40) LanC-like protein 1 LANCL1 GPR69A
Observed MW:	45 kDa
Uniprot ID :	O43813

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer using FNab04691(LANCL1 antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04691(LANCL1 antibody) at dilution of 1:300