

CD207 antibody

Product Information

Catalog No.:	FNab04692
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	CD207 molecule, langerin
Synonyms:	C-type lectin domain family 4 member K Langerin CD207 CLEC4K
Observed MW:	40 kDa
Uniprot ID :	Q9UJ71

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human skin tissue were subjected to SDS PAGE followed by western blot with FNab04692(Langerin antibody) at dilution of 1:500