

TMPO antibody

Product Information

Catalog No.:	FNab04693
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	thymopoietin
Synonyms:	Lamina-associated polypeptide 2, isoform alpha Thymopoietin isoform alpha (TP alpha) Thymopoietin-related peptide isoform alpha (TPRP isoform alpha) Thymopoietin (TP) Alternative names: Splenin Thymopentin Alternative names: TP5 TMPO LAP2
Observed MW:	51 kDa
Uniprot ID :	P42166

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin cancer using FNab04693(TMPO antibody) at dilution of 1:100



K-562 cells were subjected to SDS PAGE followed by western blot with FNab04693(TMPO antibody) at dilution of 1:500