

## LARP4 antibody

### Product Information

Catalog No.:	FNab04699
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Plays a role in the regulation of cell morphology and cytoskeletal organization.

### Immunogen information


Immunogen:	La ribonucleoprotein domain family, member 4
Synonyms:	La-related protein 4 La ribonucleoprotein domain family member 4 LARP4
Observed MW:	81 kDa
Uniprot ID :	Q71RC2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab04699(LARP4 antibody) at dilution of 1:50



Y79 cells were subjected to SDS PAGE followed by western blot with FNab04699(LARP4 antibody) at dilution of 1:300