

LCE1B antibody

Product Information

Catalog No.:	FNab04724
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Precursors of the cornified envelope of the stratum corneum.

Immunogen information


Immunogen:	late cornified envelope 1B
Synonyms:	Late cornified envelope protein 1B Late envelope protein 2 Small proline-rich-like epidermal differentiation complex protein 2A LCE1B LEP2 SPRL2A
Observed MW:	10-12 kDa
Uniprot ID :	Q5T7P3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin using FNab04724(LCE1B antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab04724(LCE1B antibody) at dilution of 1:1000