

## LCOR antibody

### Product Information

Catalog No.:	FNab04730
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May act as transcription activator that binds DNA elements with the sequence 5'-CCCTATCGATCGATCTCTACCT-3'(By similarity). Repressor of ligand-dependent transcription activation by target nuclear receptors. Repressor of ligand-dependent transcription activation by ESR1, ESR2, NR3C1, PGR, RARA, RARB, RARG, RXRA and VDR.

### Immunogen information


Immunogen:	ligand dependent nuclear receptor corepressor
Synonyms:	Ligand-dependent corepressor (LCoR) Mblk1-related protein 2 LCOR C10orf12 KIAA1795 MLR2
Observed MW:	50 kDa
Uniprot ID :	Q96JN0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab04730(LCOR antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04730(LCOR antibody) at dilution of 1:300