

LDB3 antibody

Product Information

Catalog No.:	FNab04734
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as an adapter in striated muscle to couple protein kinase C-mediated signaling via its LIM domains to the cytoskeleton.

Immunogen information


Immunogen:	LIM domain binding 3
Synonyms:	LIM domain-binding protein 3 Protein cypher Z-band alternatively spliced PDZ-motif protein LDB3 KIAA0613 ZASP
Observed MW:	32 kDa, 78-80 kDa
Uniprot ID :	O75112

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab04734(LDB3 antibody) at dilution of 1:50



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab04734(LDB3 antibody) at dilution of 1:1000