

LDHB-Specific antibody

Product Information

Catalog No.:	FNab04739
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	lactate dehydrogenase B
Synonyms:	L-lactate dehydrogenase B chain (LDH-B) LDH heart subunit (LDH-H) Renal carcinoma antigen NY-REN-46 LDHB
Observed MW:	36 kDa
Uniprot ID :	P07195

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:100-1:400; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab04739(LDHB-Specific antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04739(LDHB-Specific antibody) at dilution of 1:300