

## LDHD antibody

### Product Information

Catalog No.:	FNab04741
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	lactate dehydrogenase D
Synonyms:	Probable D-lactate dehydrogenase, mitochondrial (DLD, Lactate dehydrogenase D) LDHD
Observed MW:	45-54 kDa
Uniprot ID :	Q86WU2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human osteosarcoma using FNab04741(LDHD antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04741(LDHD antibody) at dilution of 1:100