

## LGALS3BP antibody

### Product Information

Catalog No.:	FNab04757
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Promotes integrin-mediated cell adhesion. May stimulate host defense against viruses and tumor cells.

### Immunogen information


Immunogen:	lectin, galactoside-binding, soluble, 3 binding protein
Synonyms:	Galectin-3-binding protein Basement membrane autoantigen p105 Lectin galactoside-binding soluble 3-binding protein Mac-2-binding protein (MAC2BP, Mac-2 BP) Tumor-associated antigen 90K LGALS3BP M2BP
Observed MW:	65-70 kDa
Uniprot ID :	Q08380

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-LGALS3BP (IP:FNab04757, 4ug;  
Detection:FNab04757 1:300) with HepG2 cells  
lysate 3480ug.



HepG2 cells were subjected to SDS PAGE  
followed by western blot with  
FNab04757(LGALS3BP antibody) at dilution of  
1:1000