

LGI2 antibody

Product Information

Catalog No.:	FNab04761
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	leucine-rich repeat LGI family, member 2
Synonyms:	Leucine-rich repeat LGI family member 2 LGI1-like protein 2 Leucine-rich glioma-inactivated protein 2 LGI2 KIAA1916 LGIL2
Observed MW:	62 kDa
Uniprot ID :	Q8N0V4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab04761(LGI2 antibody) at dilution of 1:300