

LGI3 antibody

Product Information

Catalog No.:	FNab04762
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May participate in the regulation of neuronal exocytosis.

Immunogen information


Immunogen:	leucine-rich repeat LGI family, member 3
Synonyms:	Leucine-rich repeat LGI family member 3 LGI1-like protein 4 Leucine-rich glioma-inactivated protein 3 LGI3 LGIL4
Observed MW:	55-62 kDa
Uniprot ID :	Q8N145

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung using FNab04762(LGI3 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04762(LGI3 antibody) at dilution of 1:1000