

EIF2D antibody

Product Information

Catalog No.:	FNab04767
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	ligatin
Synonyms:	Eukaryotic translation initiation factor 2D (eIF2d) Hepatocellular carcinoma-associated antigen 56 Ligatin EIF2D HCA56 LGTN
Observed MW:	65 kDa
Uniprot ID :	P41214

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:3000; IF: 1:10-1:100
Image:	

HeLa cells were subjected to SDS PAGE followed by western blot with FNab04767(LGTN antibody) at dilution of 1:500

