

LILRA5 antibody

Product Information

Catalog No.:	FNab04777
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in triggering innate immune responses. Does not seem to play a role for any class I MHC antigen recognition.

Immunogen information

Immunogen:	leukocyte immunoglobulin-like receptor, subfamily A(with TM domain), member 5
Synonyms:	Leukocyte immunoglobulin-like receptor subfamily A member 5 CD85 antigen-like family member F Immunoglobulin-like transcript 11 (ILT-11) Leukocyte immunoglobulin-like receptor 9 (LIR-9) LILRA5 ILT11 LILRB7 LIR9
Observed MW:	31-33 kDa
Uniprot ID :	A6NI73

Application

Reactivity:	Human
Tested Application:	ELISA, WB, FC
Recommended dilution: WB:	1:500-1:2000
Image:	



K562 cells were subjected to SDS PAGE followed by western blot with FNab04777(LILRA5 antibody) at dilution of 1:300