

PPFIA1 antibody

Product Information

Catalog No.:	FNab04797
Size:	100μg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	protein tyrosine phosphatase, receptor type, f polypeptide(PTPRF), interacting protein(liprin), alpha 1
Synonyms:	Liprin-alpha-1 LAR-interacting protein 1 (LIP-1) Protein tyrosine phosphatase receptor type f polypeptide-interacting protein alpha-1 (PTPRF-interacting protein alpha-1) PPFIA1 LIP1
Observed MW:	160 kDa
Uniprot ID :	Q13136

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab04797(PPF1A1 antibody) at dilution of 1:100



IP Result of anti-Liprin alpha 1 (IP:FNab04797, 4ug; Detection:FNab04797 1:500) with MCF-7 cells lysate 3200ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04797(PPF1A1 antibody) at dilution of 1:500