

LMF1 antibody

Product Information

Catalog No.:	FNab04805
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the maturation of specific proteins in the endoplasmic reticulum. Required for maturation and transport of active lipoprotein lipase(LPL) through the secretory pathway. Each LMF1 molecule chaperones 50 or more molecules of LPL.

Immunogen information


Immunogen:	lipase maturation factor 1
Synonyms:	Lipase maturation factor 1 Transmembrane protein 112 LMF1 C16orf26 TMEM112
Observed MW:	65-70 kDa
Uniprot ID :	Q96S06

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab04805(LMF1 Antibody) at dilution of 1:50



SKOV-3 cells were subjected to SDS PAGE followed by western blot with FNab04805(LMF1 Antibody) at dilution of 1:600