

LMF2 antibody

Product Information

Catalog No.:	FNab04806
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the maturation of specific proteins in the endoplasmic reticulum. May be required for maturation and transport of active lipoprotein lipase(LPL) through the secretory pathway(By similarity).

Immunogen information


Immunogen:	lipase maturation factor 2
Synonyms:	Lipase maturation factor 2 Transmembrane protein 112B Transmembrane protein 153 LMF2 TMEM112B TMEM153
Observed MW:	80 kDa
Uniprot ID :	Q9BU23

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney using FNab04806(LMF2 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04806(LMF2 antibody) at dilution of 1:800