

LMOD1 antibody

Product Information

Catalog No.:	FNab04808
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


The leiomodins 1 protein has a putative membrane-spanning region and 2 types of tandemly repeated blocks. The transcript is expressed in all tissues tested, with the highest levels in thyroid, eye muscle, skeletal muscle, and ovary. Increased expression of leiomodins 1 may be linked to Graves' disease and thyroid-associated ophthalmopathy.

Immunogen information

Immunogen:	leiomodins 1 (smooth muscle)
Synonyms:	Leiomodins-1 64 kDa autoantigen 1D 64 kDa autoantigen 1D3 64 kDa autoantigen D1 Leiomodins, muscle form Smooth muscle leiomodins (SM-Lmod) Thyroid-associated ophthalmopathy autoantigen LMOD1
Observed MW:	Refer to figures
Uniprot ID :	P29536

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon using FNab04808(LMOD1 antibody) at dilution of 1:100