

LMX1B antibody

Product Information

Catalog No.:	FNab04812
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Essential for the specification of dorsal limb fate at both the zeugopodal and autopodal levels.

Immunogen information

Immunogen:	LIM homeobox transcription factor 1, beta
Synonyms:	LIM homeobox transcription factor 1-beta LIM/homeobox protein 1.2 (LMX-1.2) LIM/homeobox protein LMX1B LMX1B
Observed MW:	42 kDa
Uniprot ID :	O60663

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04812(LMX1B Antibody) at dilution of 1:600