

BORCS5 antibody

Product Information

Catalog No.:	FNab04819
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	loss of heterozygosity, 12, chromosomal region 1
Synonyms:	BLOC-1-related complex subunit 5 Loss of heterozygosity 12 chromosomal region 1 Myristoylated lysosomal protein (Myrlysin) BORCS5 LOH12CR1
Observed MW:	17-22 kDa
Uniprot ID :	Q969J3

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with FNab04819(LOH12CR1 antibody) at dilution of 1:300