

LORICRIN antibody

Product Information

Catalog No.:	FNab04822
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	loricrin
Synonyms:	Loricrin LORICRIN LOR LRN
Observed MW:	26 kDa
Uniprot ID :	P23490

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:10000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human skin tissue slide using FNab04822(LOR Antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab04822(LOR antibody) at dilution of 1:500