

## ADGRL3 antibody

### Product Information

Catalog No.:	FNab04829
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	latrophilin 3
Synonyms:	Adhesion G protein-coupled receptor L3 Calcium-independent alpha-latrotoxin receptor 3 (CIRL-3) Latrophilin-3 Lectomedin-3 ADGRL3 KIAA0768 LEC3 LPHN3
Observed MW:	139 kDa
Uniprot ID :	Q9HAR2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-LPHN3 (IP:FNab04829, 4ug;  
Detection:FNab04829 1:300) with mouse brain  
tissue lysate 4000ug.



mouse brain tissue were subjected to SDS PAGE  
followed by western blot with FNab04829(LPHN3  
antibody) at dilution of 1:300