

LRP2BP antibody

Product Information

Catalog No.:	FNab04846
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May act as an adapter that regulates LRP2 function.

Immunogen information

Immunogen:	LRP2 binding protein
Synonyms:	LRP2-binding protein Megalin-binding protein (MegBP) LRP2BP KIAA1325
Observed MW:	45 kDa
Uniprot ID :	Q9P2M1

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab04846(LRP2BP antibody) at dilution of 1:300