

## LRP2-Specific antibody

### Product Information

Catalog No.:	FNab04847
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	low density lipoprotein-related protein 2
Synonyms:	Low-density lipoprotein receptor-related protein 2 (LRP-2) Glycoprotein 330 (gp330) Megalin LRP2
Observed MW:	280 kDa
Uniprot ID :	P98164

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab04847(LRP2-Specific Antibody) at dilution of 1:200



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04847(LRP2-Specific antibody) at dilution of 1:300