

LRRC10 antibody

Product Information

Catalog No.:	FNab04853
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play important roles in cardiac development and/or cardiac function.

Immunogen information

Immunogen:	leucine rich repeat containing 10
Synonyms:	Leucine-rich repeat-containing protein 10 LRRC10
Observed MW:	32 kDa
Uniprot ID :	Q5BKY1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04853(LRRC10 antibody) at dilution of 1:300