

LRRC47 antibody

Product Information

Catalog No.:	FNab04858
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	leucine rich repeat containing 47
Synonyms:	Leucine-rich repeat-containing protein 47 LRRC47 KIAA1185
Observed MW:	63 kDa
Uniprot ID :	Q8N1G4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas slide using FNab04858(LRRC47 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab04858(LRRC47 antibody) at dilution of 1:1000