

## LRRC4C antibody

### Product Information

Catalog No.:	FNab04859
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May promote neurite outgrowth of developing thalamic neurons.

### Immunogen information


Immunogen:	leucine rich repeat containing 4C
Synonyms:	Leucine-rich repeat-containing protein 4C Netrin-G1 ligand (NGL-1) LRRC4C KIAA1580 NGL1
Observed MW:	72 kDa
Uniprot ID :	Q9HCJ2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab04859(LRRC4C antibody) at dilution of 1:50



Y79 cells were subjected to SDS PAGE followed by western blot with FNab04859(LRRC4C antibody) at dilution of 1:500