

LSM10 antibody

Product Information

Catalog No.:	FNab04871
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Appears to function in the U7 snRNP complex that is involved in histone 3'-end processing. Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases. Binds specifically to U7 snRNA. Binds to the downstream cleavage product(DCP) of histone pre-mRNA in a U7 snRNP dependent manner.

Immunogen information

Immunogen:	LSM10, U7 small nuclear RNA associated
Synonyms:	U7 snRNA-associated Sm-like protein LSm10 LSM10
Observed MW:	14 kDa
Uniprot ID :	Q969L4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



human spleen tissue were subjected to SDS PAGE followed by western blot with FNab04871(LSM10 antibody) at dilution of 1:100