

LSM11 antibody

Product Information

Catalog No.:	FNab04872
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of the U7 snRNP complex that is involved in the histone 3'-end pre-mRNA processing(By similarity). Increases U7 snRNA levels but not histone 3'-end pre-mRNA processing activity, when overexpressed. Required for cell cycle progression from G1 to S phases. Binds specifically to the Sm-binding site of U7 snRNA.

Immunogen information

Immunogen:	LSM11, U7 small nuclear RNA associated
Synonyms:	U7 snRNA-associated Sm-like protein LSm11 LSM11
Observed MW:	45-50 kDa
Uniprot ID :	P83369

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04872(LSM11 antibody) at dilution of 1:600