

LUM antibody

Product Information

Catalog No.:	FNab04890
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	lumican
Synonyms:	Lumican Keratan sulfate proteoglycan lumican (KSPG lumican) LUM LDC SLRR2D
Observed MW:	65-90 kDa
Uniprot ID :	P51884

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



JurKat cells were subjected to SDS PAGE followed by western blot with FNab04890(Lumican antibody) at dilution of 1:500