

LXN antibody

Product Information

Catalog No.:	FNab04893
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Hardly reversible, non-competitive, and potent inhibitor of CPA1, CPA2 and CPA4. May play a role in inflammation.

Immunogen information


Immunogen:	latexin
Synonyms:	Latexin Endogenous carboxypeptidase inhibitor (ECI) Protein MUM Tissue carboxypeptidase inhibitor (TCI) LXN
Observed MW:	26 kDa
Uniprot ID :	Q9BS40

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-LXN (IP:FNab04893, 4ug;
Detection:FNab04893 1:500) with mouse lung
tissue lysate 5160ug.



mouse lung tissue were subjected to SDS PAGE
followed by western blot with FNab04893(LXN
Antibody) at dilution of 1:500