

LYRM7 antibody

Product Information

Catalog No.:	FNab04908
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Assembly factor required for Rieske Fe-S protein UQCRFS1 incorporation into the cytochrome b-c1(CIII) complex. Functions as a chaperone, binding to this subunit within the mitochondrial matrix and stabilizing it prior to its translocation and insertion into the late CIII dimeric intermediate within the mitochondrial inner membrane.

Immunogen information

Immunogen:	Lym7 homolog(Mouse)
Synonyms:	Complex III assembly factor LYRM7 LYR motif-containing protein 7 LYRM7 C5orf31 MZM1L
Observed MW:	None
Uniprot ID :	Q5U5X0

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, IF, FC
Recommended dilution:	IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab04908(LYRM7 antibody) at dilution of 1:50