

MAB21L2 antibody

Product Information

Catalog No.:	FNab04919
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for several aspects of embryonic development including normal development of the eye.

Immunogen information

Immunogen:	mab-21-like 2(C. elegans)
Synonyms:	Protein mab-21-like 2 MAB21L2
Observed MW:	40 kDa
Uniprot ID :	Q9Y586

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04919(MAB21L2 antibody) at dilution of 1:200