

MAD2L2 antibody

Product Information

Catalog No.:	FNab04926
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background


The protein encoded by this gene is a component of the mitotic spindle assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. The encoded protein, which is similar to MAD2L1, is capable of interacting with ADAM9, ADAM15, REV1, and REV3 proteins.

Immunogen information


Immunogen:	MAD2 mitotic arrest deficient-like 2 (yeast)
Synonyms:	Mitotic spindle assembly checkpoint protein MAD2B Mitotic arrest deficient 2-like protein 2 (MAD2-like protein 2) REV7 homolog (hREV7) MAD2L2 MAD2B REV7
Observed MW:	24 kDa
Uniprot ID :	Q9UI95

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human gastric cancer using FNab04926(MAD2L2 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04926(MAD2L2 antibody) at dilution of 1:1000