

MAGEA10 antibody

Product Information

Catalog No.:	FNab04935
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Not known, though may play a role in embryonal development and tumor transformation or aspects of tumor progression.

Immunogen information

Immunogen:	melanoma antigen family A, 10
Synonyms:	Melanoma-associated antigen 10 Cancer/testis antigen 1.10 (CT1.10) MAGE-10 antigen MAGEA10 MAGE10
Observed MW:	41 kDa, 62 kDa
Uniprot ID :	P43363

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04935(MAGEA10 antibody) at dilution of 1:500