

MAGEA11 antibody

Product Information

Catalog No.:	FNab04936
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Acts as androgen receptor coregulator that increases androgen receptor activity by modulating the receptors interdomain interaction. May play a role in embryonal development and tumor transformation or aspects of tumor progression.

Immunogen information

Immunogen:	melanoma antigen family A, 11
Synonyms:	Melanoma-associated antigen 11 Cancer/testis antigen 1.11 (CT1.11) MAGE-11 antigen MAGEA11 MAGE11
Observed MW:	48 kDa, 44 kDa
Uniprot ID :	P43364

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



A375 cells were subjected to SDS PAGE followed by western blot with FNab04936(MAGEA11 antibody) at dilution of 1:200