

MAGEB1 antibody

Product Information

Catalog No.:	FNab04941
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	melanoma antigen family B, 1
Synonyms:	Melanoma-associated antigen B1 Cancer/testis antigen 3.1 (CT3.1) DSS-AHC critical interval MAGE superfamily 10 (DAM10) MAGE-B1 antigen MAGE-XP antigen MAGEB1 MAGEL1 MAGEXP
Observed MW:	43-45 kDa
Uniprot ID :	P43366

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04941(MAGEB1 antibody) at dilution of 1:2000