

MAGEC3 antibody

Product Information

Catalog No.:	FNab04945
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	melanoma antigen family C, 3
Synonyms:	Melanoma-associated antigen C3 Cancer/testis antigen 7.2 (CT7.2) Hepatocellular carcinoma-associated antigen 2 MAGE-C3 antigen MAGEC3 HCA2
Observed MW:	70-72 kDa
Uniprot ID :	Q8TD91

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer using FNab04945(MAGEC3 antibody) at dilution of 1:50



A375 cells were subjected to SDS PAGE followed by western blot with FNab04945(MAGEC3 Antibody) at dilution of 1:600