

## MAGED2 antibody

### Product Information

Catalog No.:	FNab04948
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	melanoma antigen family D, 2
Synonyms:	Melanoma-associated antigen D2 11B6 Breast cancer-associated gene 1 protein (BCG-1) Hepatocellular carcinoma-associated protein JCL-1 MAGE-D2 antigen MAGED2 BCG1
Observed MW:	65 kDa
Uniprot ID :	Q9UNF1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IF: 1:20-1:200
Image:	



IP Result of anti-MAGED2 (IP:FNab04948, 4ug;  
Detection:FNab04948 1:1000) with HepG2 cells  
lysate 3200ug.



A431 cells were subjected to SDS PAGE followed  
by western blot with FNab04948(MAGED2  
antibody) at dilution of 1:1000