

## MAGED4 antibody

### Product Information

Catalog No.:	FNab04949
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	melanoma antigen family D, 4
Synonyms:	Melanoma-associated antigen D4 MAGE-D4 antigen MAGE-E1 antigen MAGED4 KIAA1859 MAGED4A MAGEE1
Observed MW:	42-60 kDa
Uniprot ID :	Q96JG8

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:1000-1:5000
Image:	



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab04949(MAGED4 antibody) at dilution of 1:5000