

MAK16 antibody

Product Information

Catalog No.:	FNab04957
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	MAK16 homolog
Synonyms:	Protein MAK16 homolog NNP78 Protein RBM13 MAK16 RBM13
Observed MW:	48 kDa
Uniprot ID :	Q9BXY0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:1000-1:4000; IF: 1:50-1:500
Image:	



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using FNab04957 (MAK16 antibody) at dilution of 1:100 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04957(MAK16 antibody) at dilution of 1:500