

## MAML3 antibody

### Product Information

Catalog No.:	FNab04961
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Acts as a transcriptional coactivator for NOTCH proteins. Has been shown to amplify NOTCH-induced transcription of HES1.

### Immunogen information

Immunogen:	mastermind-like 3(Drosophila)
Synonyms:	Mastermind-like protein 3 (Mam-3) MAML3 KIAA1816
Observed MW:	122 kDa
Uniprot ID :	Q96JK9

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04961(MAML3 antibody) at dilution of 1:2000