

## SCGB2A2 antibody

### Product Information

Catalog No.:	FNab04964
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4D8
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	secretoglobin, family 2A, member 2
Synonyms:	Mammaglobin-A Mammaglobin-1 Secretoglobin family 2A member 2 SCGB2A2 MGB1 UGB2
Observed MW:	6-10 kDa
Uniprot ID :	Q13296

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer slide using FNab04964(SCGB2A2 Antibody) at dilution of 1:50