

MAN1A2 antibody

Product Information

Catalog No.:	FNab04965
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the maturation of Asn-linked oligosaccharides. Progressively trim alpha-1,2-linked mannose residues from Man(9)GlcNAc(2) to produce Man(5)GlcNAc(2).

Immunogen information


Immunogen:	mannosidase, alpha, class 1A, member 2
Synonyms:	Mannosyl-oligosaccharide 1,2-alpha-mannosidase IB Mannosidase alpha class 1A member 2 Processing alpha-1,2-mannosidase IB (Alpha-1,2-mannosidase IB) MAN1A2 MAN1B
Observed MW:	73 kDa
Uniprot ID :	O60476

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta using FNab04965(MAN1A2 antibody) at dilution of 1:50



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab04965(MAN1A2 antibody) at dilution of 1:200