

MANBAL antibody

Product Information

Catalog No.:	FNab04967
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mannosidase, beta A, lysosomal-like
Synonyms:	Protein MANBAL MANBAL
Observed MW:	10 kDa
Uniprot ID :	Q9NQG1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500 - 1:5000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab04967(MANBAL antibody) at dilution of 1:100



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04967(MANBAL Antibody) at dilution of 1:1000