

MAP3K4 antibody

Product Information

Catalog No.:	FNab04978
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Component of a protein kinase signal transduction cascade. Activates the CSBP2, P38 and JNK MAPK pathways, but not the ERK pathway. Specifically phosphorylates and activates MAP2K4 and MAP2K6.

Immunogen information

Immunogen:	mitogen-activated protein kinase kinase kinase 4
Synonyms:	Mitogen-activated protein kinase kinase kinase 4 MAP three kinase 1 MAPK/ERK kinase kinase 4 (MEK kinase 4, MEKK 4) MAP3K4 KIAA0213 MAPKKK4 MEKK4 MTK1
Observed MW:	170-180 kDa
Uniprot ID :	Q9Y6R4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



K-562 cells were subjected to SDS PAGE followed by western blot with FNab04978(MAP3K4 antibody) at dilution of 1:2000